



## Product Datasheet

### Human IGFBP-4 Protein (orb348805)

<b>Description</b>	Insulin-like growth factor-binding protein 4 (IGFBP-4) is also known as IGF-binding protein 4, which...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	IGFBP-4
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 27.9 kDa. The protein migrates as 33-37 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Purity</b>	95%
<b>Source</b>	Human IGFBP-4, His Tag (orb348805) is expressed from human 293 cells (HEK293). It contains AA Asp 22 - Glu 258 (Accession # P22692-1).
<b>MW</b>	27.9 kDa
<b>NCBI</b>	<a href="#">NP_001543.2</a>
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human IGFBP-4 prote...



Immobilized Human IGFBP-4, His Tag (Cat....