


Product Datasheet

Zika virus NS1 Protein (orb348810)

Description	Zika virus (ZIKV) is a member of the virus family Flaviviridae and the genus Flavivirus, transmitted...
Reactivity	Virus
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	NS1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 42.0 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Val 796 - Ser 1148
Purity	95%
Source	Zika virus NS1, His Tag (orb348810) is expressed from human 293 cells (HEK293). It contains AA Val 796 - Ser 1148 (Accession # ALU33341).
MW	42.0 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Zika virus
NS1 prot...