

Product Datasheet

C20orf197 antibody (orb34882)

Description	C20orf197 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Uncharacterized protein C20orf197 protein (1-126AA)
Target	C20orf197
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N268
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.