


Product Datasheet

Human Galectin-9 / LGALS9 Protein (orb348820)

Description	Galectin-9 (LGALS9) is also known as tumor antigen HOM-HD-21 and ecalectin. Galectins are a family o...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Galectin-9
Preservatives	20 mM MOPS, 50 mM NaCl, 1 mM EDTA, 2 mM DTT, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	N-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 36.6 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	BAA31542.1
Purity	95%
Source	Human Galectin-9, His Tag (orb348820) is expressed from human 293 cells (HEK293). It contains AA Ala 2 - Thr 323 (Accession # O00182-2).
MW	36.6 kDa
NCBI	BAA31542.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Galectin-9 pr...