



Product Datasheet

CD52 antibody (orb348863)


Description	Human monoclonal antibody to CD52
Species/Host	Human
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, WB
Immunogen	Human lymphocytes.
Target	CD52
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P31358
Expiration Date	12 months from date of receipt.



IF image of rat thymus tissue using CD52...



Immunofluorescence analysis of rat thymu...



Immunohistochemical staining of rat sple...