



Product Datasheet

TNF alpha antibody (orb348870)


Description	Human monoclonal antibody to TNF alpha
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human TNF-alpha
Target	TNF alpha
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	The antibody binds specifically to TNF-alpha, a pro-inflammatory cytokine which stimulates release of other cytokines e.g. IL-1, IL-6 and GM-CSF as well as secretion of ICAM-1 and VCAM-1 from endothelial cells. It stimulates acute phase responses and acts as a chemoattractant for NK cells and monocytes/macrophages. The antibody reached Phase II clinical trials with Celltech, UCB.
Isotype	IgG4
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P01375
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human th...



Immunohistochemical staining of paraffin...



IHC-P staining of rat spinal cord tissue...