

Product Datasheet

NgR2 antibody (orb348888)

Description	Rabbit monoclonal antibody to NgR2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IP, WB
Immunogen	(His)6-tagged fragment of rat NgR2 (amino acids 279-420).
Target	NgR2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	The antibody binds specifically to residues 351-399 of NgR2, a member of the Nogo-66 receptor family which inhibits axonal growth and sporulating in the CNS.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q86UN3
Expiration Date	12 months from date of receipt.