

Product Datasheet

CD38 antibody (orb348893)

Description	Rabbit monoclonal antibody to CD38
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Viable cells (1 x 10E7) of the human cell line, BJAB (Human Burkitt lymphoma B cell line), which were injected intraperitoneally into 8-week-old BALB/c mice.
Target	CD38
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	A customer reports successfully using this antibody (coupled to magnetic beads) for isolation of plasma cells from the intestinal lamina propria.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P28907
Expiration Date	12 months from date of receipt.