

Product Datasheet

P-glycoprotein antibody (orb348896)

Description	Rabbit monoclonal antibody to P-glycoprotein
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IP
Immunogen	Drug resistant K-562 cells.
Target	P-glycoprotein
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P08183
Expiration Date	12 months from date of receipt.