

Product Datasheet

CD3 epsilon antibody (orb348901)

Description	Rabbit monoclonal antibody to CD3E
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	H-2Kb-specific murine cytotoxic T-lymphocyte (CTL) clone.
Target	CD3 epsilon
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P22646
Expiration Date	12 months from date of receipt.