

## Product Datasheet

### IL-2R antibody (orb348909)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit monoclonal antibody to IL-2R                                              |
| <b>Species/Host</b>        | Human                                                                            |
| <b>Reactivity</b>          | Human, Primate                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                     |
| <b>Tested Applications</b> | Blocking, FC                                                                     |
| <b>Immunogen</b>           | Cultured T-cells derived from a patient with Mycosis fungoides.                  |
| <b>Target</b>              | IL-2R                                                                            |
| <b>Preservatives</b>       | PBS with 0.02% Proclin 300.                                                      |
| <b>Concentration</b>       | 1 mg/ml                                                                          |
| <b>Storage</b>             | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                       |
| <b>Purity</b>              | Purified                                                                         |
| <b>Uniprot ID</b>          | <b>P01589</b>                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |