

Product Datasheet

Integrin beta 1 antibody (orb348910)

Description	Rabbit monoclonal antibody to Integrin beta 1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Immunogen	Human T-cell lymphoma cells injected into Biozzi high-responder mice.
Target	CD29
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to CD29.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P05556
Expiration Date	12 months from date of receipt.