



Product Datasheet

Smad2 (phospho-Ser118) antibody (orb34892)

Description	Rabbit polyclonal antibody against Smad2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IF
Immunogen	Synthetic Peptide
Target	SMAD2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52306
Uniprot ID	Q15796
NCBI	NP_005892.1 , NP_001003652.1
Dilution Range	IF/ICC: 1:10-50, DOT: 1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...