

Product Datasheet

CD11a antibody (orb348923)

Description	Rabbit monoclonal antibody to CD11a
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	Efalizumab binds to CD11a, a subunit of leukocyte function antigen-1 (LFA-1), which is expressed on all leukocytes.
Target	CD11a
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	Efalizumab binds to the LFA-1 alpha I domain blocking the ICAM-1 binding site via steric hindrance.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P20701
Expiration Date	12 months from date of receipt.