


Product Datasheet

Human Caspase 8/CASP8 ELISA Kit (orb348939)

Description	Human Caspase 8/CASP8 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S217-D384 & L385-D479
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	55391 MW
Uniprot ID	Q14790
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is around 120 pg/ml.
Expiration Date	Please enquire.



Human
Caspase 8
PicoKine
ELISA Kit
stand...