

## Product Datasheet

# Human AGRP/ART/Agouti-related protein ELISA Kit (orb348940)

**Description** Human AGRP/ART/Agouti-related protein ELISA Kit

**Reactivity** Human

**Range** 7.8 pg/ml - 500 pg/ml

**Tested Applications** ELISA

**Immunogen** Expression system for standard: NS0; Immunogen sequence: A21-T132

**Storage** Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

**Note** For research use only


**Sample Types** cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.

**MW** 14440 MW

**Uniprot ID** **000253**

**Sensitivity** 4 pg/ml

**Expiration Date** Please enquire.



Human  
AGRP  
PicoKine  
ELISA Kit  
standard  
C...