


Product Datasheet

Human Relaxin 2 ELISA Kit (orb348941)

Description	Human Relaxin 2 ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q162-C185 (A) & D25-S53 (B)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	21043 MW
Uniprot ID	P04090
Sensitivity	4 pg/ml
Expiration Date	Please enquire.



Human
Relaxin 2
PicoKine
ELISA Kit
stand...