


Product Datasheet

Rat Eotaxin ELISA Kit (orb348946)

Description	Rat Eotaxin ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: H24-P97
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (EDTA, citrate).
MW	10851 MW
Uniprot ID	P97545
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:4, concentration in serum and plasma is around 1 ng/ml.
Expiration Date	Please enquire.



Rat Eotaxin
PicoKine
ELISA Kit
standard ...