


Product Datasheet

Human APOE/Apolipoprotein E ELISA Kit (orb348956)

Description	Human APOE/Apolipoprotein E ELISA Kit
Reactivity	Human
Range	3.12 ng/ml - 200 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K19-H317
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, milk and urine.
MW	36154 MW
Uniprot ID	P02649
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1000, concentration in serum and plasma is around 20-40 ug/ml.
Expiration Date	Please enquire.



Human
APOE
PicoKine
ELISA Kit
standard
C...