


Product Datasheet

Human APOA1/Apolipoprotein A-I ELISA Kit (orb348957)

Description	Human APOA1/Apolipoprotein A-I ELISA Kit
Reactivity	Human
Range	3.12 ng/ml - 200 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: D25-Q267
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	30778 MW
Uniprot ID	P02647
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:200000, concentration in serum and plasma is around 5mg/ml.
Expiration Date	Please enquire.



Human
APOA1
PicoKine
ELISA Kit
standard ...