

Product Datasheet

IL10 antibody (orb348958)

Description

IL10 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human IL10. Please contact us for the exact immunogen sequence. The peptide is available as orb374965.

Target

IL10

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

21 kDa

Uniprot ID

[P22301](#), [P29456](#), [P18893](#)

NCBI

[21](#), [000572](#), [000563](#), [11](#)

Entrez



[3586](#)

Dilution Range


IHC-P: 1:100-500, WB: 1:100-2000, (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



IHC-P image of rat lymph node tissue
usi...



Western blot analysis
of rat lymph node ...



Immunohistochemical
staining of rat colo...