

Product Datasheet

VCAM1 antibody (orb348962)

Description

VCAM1 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human VCAM1. Please contact us for the exact immunogen sequence. The peptide is available as orb374969.

Target

VCAM1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

81 kDa

Uniprot ID

[P29534](#), [P19320](#), [P29533](#)

NCBI

[001069](#), [001078](#)

Entrez


[7412](#)

Dilution Range


IHC-P: 1:200, WB: 1:100-2000, (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of Mouse lu...



Immunohistochemical staining of paraffin...



WB analysis of rat testis (lane 1), rat ...