

Product Datasheet

CD3 antibody (orb348965)

Description

CD3 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human CD3. Please contact us for the exact immunogen sequence. The peptide is available as orb374972.

Target

CD3E

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

19 kDa

Uniprot ID

[P22646](#), [P07766](#)

NCBI


[000733](#), [31](#), [000724](#), [11](#)

Entrez
[916](#)
Dilution Range


IHC-P: 1:100-1000, WB: 1:100-1000, (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of rat bone...



IHC-P image of rat skin tissue using CD3...