




Product Datasheet

ASCL1 antibody (orb348967)

Description	Rabbit polyclonal antibody to ASCL1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	ASCL1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25454
Uniprot ID	P50553
NCBI	NP_004307.2
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis of 293 cell lysates using AS...



WB analysis of human brain lysate (Lane1...