

Product Datasheet

PPM1K antibody (FITC) (orb34897)

Description	PPM1K antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Protein phosphatase 1K, mitochondrial protein (30-372AA)
Target	PPM1K
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N3J5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.