




Product Datasheet

BRAF antibody (orb348970)

Description	Mouse monoclonal antibody to BRAF.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant Protein
Target	BRAF
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1κ
Clonality	Monoclonal
MW	84437
Uniprot ID	P15056
NCBI	NP_004324.2
Dilution Range	WB: 1:500-16000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



WB analysis of MCF-7 whole cell lysate (...)



Immunofluorescence analysis of C2C12 cells...