



Product Datasheet

CD71 antibody (orb348973)

Description	Mouse monoclonal antibody to CD71.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TFRC
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	84871
Uniprot ID	P02786
NCBI	NP_003225.2 , NP_001121620.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.



WB analysis of HT-1080 whole cell lysate...