

Product Datasheet

Syndecan 1 antibody (orb348974)

Description	Mouse monoclonal antibody to Syndecan 1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	SDC1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	32462
Uniprot ID	P18827
NCBI	NP_002988.3 , NP_001006947.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

