



Product Datasheet

POTE antibody (orb348975)

Description	Rabbit polyclonal antibody to POTE.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	POTED
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	66394
Uniprot ID	Q86YR6
NCBI	NP_001244291.1 , NP_778146.2
Dilution Range	WB: 1:500, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pr...