




Product Datasheet

SIGLEC11 antibody (orb348978)

Description	Rabbit polyclonal antibody to SIGLEC11.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	SIGLEC11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	75795
Uniprot ID	Q96RL6
NCBI	NP_001128635.1 , NP_443116.2
Dilution Range	WB: 1:500, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



WB analysis of
HepG2 whole cell
lysate u...



Immunohistochemical
staining of human
lu...