




Product Datasheet

MYH9 antibody (orb348979)

Description	Rabbit polyclonal antibody to MYH9.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	MYH9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	226532
Uniprot ID	P35579
NCBI	NP_002464.1
Dilution Range	WB: 1:1000, WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.



WB analysis
of A431
whole cell
lysate (L...



Flow
cytometric
analysis of
HeLa cells u...