



Product Datasheet

RGS19 antibody (orb348982)

Description	Rabbit polyclonal antibody to RGS19.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	RGS19
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	24636
Uniprot ID	P49795
NCBI	NP_001034556.1 , NP_005864.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lung tissue (H-lung) showing RGS19 expression.