



Product Datasheet

Bad antibody (orb348984)

Description	Rabbit polyclonal antibody to Bad.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Bad
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22080
Uniprot ID	Q61337
NCBI	NP_031548.1
Dilution Range	WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis of A20 whole cell lysate (La...