



Product Datasheet

AGL antibody (orb348999)

Description	Rabbit polyclonal antibody to AGL.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	AGL
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	174764
Uniprot ID	P35573
NCBI	NP_000019.2, NP_000633.2, NP_000634.2, NP_000635.2, NP_000637.2, NP_000636.2
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...