

Product Datasheet

N antibody (orb349001)

Description	Rabbit polyclonal antibody to N.
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	Synthetic Peptide
Target	N
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46025
Uniprot ID	P59595
NCBI	NP_828858.1
Expiration Date	12 months from date of receipt.