




Product Datasheet

VEGFA antibody (orb349003)

Description	Rabbit polyclonal antibody to VEGFA.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	VEGFA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27042
Uniprot ID	P15692
NCBI	NP_001165097.1, NP_001020541.2, NP_001165099.1, NP_001020537.2, NP_001020540.2, NP_001165093.1, NP_001165094.1, NP_001165096.1, NP_001020539.2, NP_001165095.1, NP_001028928.1, NP_001165098.1, NP_001020538.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



WB analysis of MCF-7 whole cell lysate (...)



WB analysis of MCF-7 whole cell lysate (...)