



Product Datasheet

APOC3 antibody (orb349008)

Description	Rabbit polyclonal antibody to APOC3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	APOC3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	10852
Uniprot ID	P02656
NCBI	NP_000031.1
Dilution Range	WB: 1:2000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human he...