




Product Datasheet

PGM1 antibody (orb349024)

Description	Rabbit polyclonal antibody to PGM1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	PGM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	61449
Uniprot ID	P36871
Dilution Range	WB: 1:2000, IHC-P: 1:25
Expiration Date	12 months from date of receipt.



WB analysis
of mouse
heart lysates
(Lane...)



WB analysis
of mouse
brain lysates
(Lane...)