



Product Datasheet

PPP1R1B antibody (orb349028)

Description	Rabbit polyclonal antibody to PPP1R1B.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	PPP1R1B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22963
Uniprot ID	Q9UD71
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis
of mouse
brain lysate
using ...