



Product Datasheet

AKR7L antibody (orb349029)

Description	Rabbit polyclonal antibody to AKR7L.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	AKR7L
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36970
Uniprot ID	Q8NHP1
Dilution Range	IF: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cells...