



Product Datasheet

SRY antibody (orb349035)

Description	Rabbit polyclonal antibody to SRY.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	SRY
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23884
Uniprot ID	Q05066
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis
of HeLa whole
cell lysates
(...)