




Product Datasheet

CTCF antibody (orb349037)

Description	Rabbit polyclonal antibody to CTCF.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	CTCF
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	82785
Uniprot ID	P49711
Dilution Range	WB: 1:2000, WB: 1:2000, IHC-P: 1:25
Expiration Date	12 months from date of receipt.



WB analysis of MCF-7 whole cell lysates ...



Immunohistochemical staining of human CO...