



Product Datasheet

DPYSL5 antibody (orb349097)

Description	Mouse monoclonal antibody to DPYSL5.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	DPYSL5
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	61421
Uniprot ID	Q9BPU6
Dilution Range	IF: 1:25, WB: 1:4000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SH-SY5Y cells.