



Product Datasheet

HCK antibody (orb349100)

Description	Mouse monoclonal antibody to HCK.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	HCK
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	59600
Uniprot ID	P08631
Dilution Range	IF/ICC: 1:25, WB: 1:2000-1:4000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...