

Product Datasheet

NGEF antibody (orb349146)

Description	Rabbit polyclonal antibody to NGEF.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	NGEF
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	82496
Uniprot ID	Q8N5V2
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.

