



Product Datasheet

BMI1 antibody (orb349181)

Description	Rabbit polyclonal antibody to BMI1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	BMI1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36949
Uniprot ID	P35226
Dilution Range	IF/ICC: 1:25, WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Human Cer...