

Product Datasheet

KIF1B antibody (orb349205)

Description	Mouse monoclonal antibody to KIF1B.
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	This KIF1B antibody is generated from a mouse immunized with a recombinant protein.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
Uniprot ID	O60333
Dilution Range	WB: 1:2000, WB: 1:2000
Expiration Date	12 months from date of receipt.

