

Product Datasheet

LINGO (phospho-Ser596) antibody (orb34921)

Description	Rabbit polyclonal to LINGO (phospho-Ser596)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	Lingo1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	69101
Uniprot ID	Q9D1T0
NCBI	NP_851419.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

