



Product Datasheet

NGFR antibody (orb349211)

Description	Rabbit polyclonal antibody to NGFR.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	NGFR
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	45183
Uniprot ID	P08138
Dilution Range	FACS: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
SH-SY5Y
cell...