




Product Datasheet

LCK antibody (orb349212)

Description	Mouse monoclonal antibody to LCK.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	LCK
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	58001
Uniprot ID	P06239
Dilution Range	IF/ICC: 1:25, FACS: 1:25, IHC-P: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human th...



Flow cytometric analysis of Jurkat cells...